



PROACTIVE SOLUTIONS FOR MOLD

If you remove the moisture you can control mold problems. Here are products to help:

Ebac BD 150CR Dehumidifier

The Ebac BD 150 Building Dryer and Dehumidifier is a rugged, yet mobile piece of equipment that will operate under extreme conditions and pull large amounts of moisture from the air. In high humidity cases, the Ebac BD 150 has removed almost 27 gallons of water in a 24-hour period! The Ebac BD 150 will dehumidify a maximum of 14,300 cubic feet.



Ebac Triton Dehumidifier

The Ebac Triton Professional Dehumidifier Is a compact rugged unit. It has the capacity to handle several room-size areas at a time (13 gallons of water in a 24-hour period). Standard 115 v. plugs in anywhere. Features include rugged wheels and integral handle for mobility, automatic pump-out (20' vertical lift, 100' + horizontal), pump purge switch, heavy duty steel chassis, powder coated for durability, Ebac's "HOT GAS" defrosting feature which automatically eliminates frost build-up providing effective operation at low ambient temperatures and Ebac's new 5 year warranty.



FIBERMASTER SANITIZING SHAMPOO

Contains "twin chain quaternary ammonium compounds to reduce odor forming bacteria in carpeting. Use with portable and truck mount extractors and rotary shampoos. Use on washable synthetics. Clear color, characteristic odor.



Fresh Mint Disinfectant
Broad spectrum bactericidal activity. 2% active quaternary. Disinfecting (780 ppm): 4 oz/ gal.

TO REDUCE THE POTENTIAL FOR MOLD GROWTH, ALWAYS CLEAN CARPETS ON LOW HUMIDITY DAYS AND USE DEHUMIDIFIERS 24 HOURS BEFORE AND AFTER.



Nobles Cyclone Dryer
Keeps air moving to accelerate drying



Vaporshark
eliminates odors in areas up to 50,000 cubic feet

Dry foam and powder cleaning systems reduce moisture while deep cleaning



HOST Liberator ExtractorVac

provides deep carpet cleaning in large areas with greater dry soil removal than conventional vacuums and superior removal of HOST dry carpet cleaner.



Von Shrader LMX

Dry foam self-contained automatic extractor

CONTACT US

MD Stetson Co., Inc.
92 York Avenue
Randolph, MA
02368

www.mdstetson.com

781-986-6161
800-255-8651
Fax: 781-961-1764