

CALTECH  
INDUSTRIES, INC.

2420 Schuette Drive  
Midland, MI 48642

## DISPATCH® Hospital Cleaner Disinfectant with Bleach

### Rapid-Acting, One-Step Cleaner Disinfectant

DISPATCH® Hospital Cleaner Disinfectant with Bleach is a rapid-acting cleaner disinfectant that works in one step. A unique, stabilized bleach and detergent solution (bleach solution strength similar to the 1:10 solution strength recommended by the CDC for effective disinfecting), DISPATCH® offers bleach strength disinfecting without the time-consuming task of mixing bleach; it also offers superior cleaning capability.

### Fast-Acting Efficacy, EPA-Registered

This product has demonstrated efficacy for a broad range of organisms in one minute. The results of efficacy tests, conducted in accordance with the test methods of the Association of Official Analytical Chemists in the presence of organic soil (5% blood serum), are described below:

<b>Organism</b>	<b>Contact Time</b>
<i>Mycobacterium tuberculosis</i> (TB) (Two minutes for DISPATCH® Towels)	30 seconds
<i>Staphylococcus aureus</i>	1 minute
<i>Salmonella choleraesuis</i>	1 minute
<i>Streptococcus pyogenes</i>	1 minute
<i>Pseudomonas aeruginosa</i>	1 minute
Methicillin resistant <i>Staphylococcus aureus</i> (MRSA) (DISPATCH® Liquid only)	1 minute
Vancomycin resistant <i>Enterococcus faecalis</i> (VRE) (DISPATCH® Liquid only)	1 minute
HIV-1 (AIDS)	1 minute
Poliovirus (Type-1, Mahoney)	1 minute
Herpes Simplex (HSV-2)	1 minute
Rotavirus (DISPATCH® Liquid only)	1 minute
<i>Trichophyton mentagrophytes</i>	1 minute

### Meets Bloodborne Pathogen Guidelines

This product meets OSHA, CDC, NCCLS and Universal Precautions Guidelines for healthcare worker protection from bloodborne pathogen contamination.

### Cleans & Disinfects in One Step

Excellent for removing blood and blood products, body fluids and excretions and everyday soil  
Cleans and disinfects in one step; no precleaning required; rarely needs rinsing

### Convenient

Available in liquid (sprays or refills) and pre-moistened, disposable towels

### Versatile, Multiple Use Product

Especially useful for laboratory countertops, point-of-care testing equipment such as blood glucose monitors, IV poles, respiratory and physical therapy equipment as well as general surfaces such as carts, bed frames and springs

Excellent for cleaning and disinfecting laminated plastic countertops, stainless steel, plastic, ceramic tile, grout, vinyl plastics, fiberglass and porcelain

### **Outperforms Conventional Bleach Solutions**

Unlike hand-diluted bleach solutions, this product:

Is convenient and ready-to-use; no mixing or discarding of unused portions

Does not rapidly deteriorate; stays stable to expiration date

Is economical; eliminates daily mixing and use of a separate cleaning product

Contains a pleasant fragrance and anti-corrosive agent

**EPA Registration No. 56392-7**

**EPA Registration No. 56392-8**

CALTECH  
INDUSTRIES, INC.

2420 Schuette Drive  
Midland, MI 48642